



Independent Bariatric  
Equipment Rental Specialists  
Tel 0800 193 2523  
info@fivemobility.co.uk  
www.fivemobility.co.uk

## Bariatric Pro Plus Dynamic Mattress



The Apollo Pro Plus Bariatric dynamic mattress benefits from a uniquely designed cell configuration in order to ensure the best possible pressure relieving qualities. The cell technology offers a revolutionary 'Very High Risk' dynamic mattress, ensuring optimum patient comfort and safety.

### Rental Details;

Product Code	5MBARIPPM
SWL	318.8kg / 50 Stone
Risk Factor	Very High Risk
Cell Configuration	3 Static, 10 figure of 8 cells
Operating Cycle	Option of either 10,15,20,25 minutes
Pressure Range	20-60 mmHg
Mattress Dimensions (inflated)	200x120x20cm
Mattress Weight	8kg
Power Unit Dimensions	37x13.5x25cm
Power Unit Weight	4.2kg

## **Key Features**

- 50 stone weight limit
- Very high risk replacement mattress
- 'Cell in Cell' replacement system with static base
- Adjustable weight settings
- Timed static switch to assist in patient transfer and for good nursing practise
- Auto-lockout function
- Designed for users at very high risk, or have existing tissue damage (pressure ulcer grade 1-4)
- Manual adjustment for to user's body shape and position for optimum pressure care
- Dynamic therapy for alternating action for lower interface pressures and gentle capillary stimulation
- Humidity management to help prevention of tissue breakdown and skin maceration
- Enhanced shear and friction protection
- Autofirm option to inflate the mattress to maximum pressure to allow safe patient transfer
- Alternating cycle time variances to offer different periods of pressure relief
- Audio and visual alarms to differentiate between low pressure, power failure and system failure errors
- Seat switch to increase the pressure in the sacrum area for when the patient is sat upright
- Welded vapour permeable multi stretch PU cover to prevent fluid ingress and internal mattress contamination
- Clinically effective against bacteria / viruses with proven resistance to MRSA
- Fire retardant to Standards BS7175 Crib 5
- CPR function to ensure rapid deflation of mattress
- The mattress cover can be wipe down or completely removed and laundered up to 95°C to aid decontamination protocols

**Note: This Product is Also Available for Purchase**